

EDF-4053

## Method 23 Internal Standard Stock Solution

1.2 mL  
in Nonane

Labeled Dioxins & Furans	Concentration (ng/mL)
2,3,7,8-TCDD ( <sup>13</sup> C <sub>12</sub> ,99%)	1000
1,2,3,7,8-PeCDD ( <sup>13</sup> C <sub>12</sub> ,99%)	1000
1,2,3,6,7,8-HxCDD ( <sup>13</sup> C <sub>12</sub> ,99%)	1000
1,2,3,4,6,7,8-HpCDD ( <sup>13</sup> C <sub>12</sub> ,99%)	1000
OCDD ( <sup>13</sup> C <sub>12</sub> ,99%)	2000
2,3,7,8-TCDF ( <sup>13</sup> C <sub>12</sub> ,99%)	1000
1,2,3,7,8-PeCDF ( <sup>13</sup> C <sub>12</sub> ,99%)	1000
1,2,3,6,7,8-HxCDF ( <sup>13</sup> C <sub>12</sub> ,99%)	1000
1,2,3,4,6,7,8-HpCDF ( <sup>13</sup> C <sub>12</sub> ,99%)	1000