

EC-5453 Rapid PCB Screening Native PAR Solution
(unlabeled)

0.5 mL
in Isooctane

Unlabeled Compounds	IUPAC#	Concentration (ng/mL)
2,4,4'-TriCB	28	1000
2,2',5'-TriCB	18	1000
2,2',3,5'-TetraCB	44	1000
2,3',4',5'-TetraCB	70	1000
2,2',5,5'-TetraCB	52	1000
2,2',4,5,5'-PentaCB	101	1000
2,3,3',4',6'-PentaCB	110	1000
2,3',4,4',5'-PentaCB	118	1000
2,2',3,4',5',6'-HexaCB	149	1000
2,2',3,4,4',5'-HexaCB	138	1000
2,2',4,4',5,5'-HexaCB	153	1000
2,2',3,4,4',5,5'-HeptaCB	180	1000
2,2',3,4',5,5',6'-HeptaCB	187	1000