



Cambridge Isotope Laboratories, Inc.

3 Highwood Drive, Tewksbury, MA 01876 USA

PH: 800.322.1174 (N. America) PH: 978-749-8000

Fax: 978-749-2768 WEB: www.isotope.com

EC-5614

EPA Method 1628 Native Standard Solution

5 mL in acetonitrile

Unlabeled	IUPAC	Concentration (ng/mL)
2-MonoCB	1	80
4-MonoCB	3	80
2,2'-DiCB	4	80
2,4'-DiCB	8	80
3,3'-DiCB	11	80
4,4'-DiCB	15	80
2,2',5-TriCB	18	80
2,2',6-TriCB	19	80
2,4,4'-TriCB	28	80
2,4',5-TriCB	31	80
3,4,4'-TriCB	37	80
2,2',3,5'-TetraCB	44	80
2,2',5,5'-TetraCB	52	80
2,2',6,6'-TetraCB	54	80
2,3,4',6-TetraCB	64	80
2,3',4,4'-TetraCB	66	80
2,3',4',5-TetraCB	70	80
2,4,4',5-TetraCB	74	80
3,3',4,4'-TetraCB	77	80
2,2',3,4,4'-PentaCB	85	80
2,2',3,5',6-PentaCB	95	80
2,2',4,4',5-PentaCB	99	80
2,2',4,5,5'-PentaCB	101	80
2,2',4,6,6'-PentaCB	104	80
2,3,3',4,4'-PentaCB	105	80
2,3,3',4',6-PentaCB	110	80
2,3',4,4',5-PentaCB	118	80
3,3',4,4',5-PentaCB	126	80

Unlabeled	IUPAC	Concentration (ng/mL)
2,2',3,3',4,6'-HexaCB	132	80
2,2',3,4,4',5'-HexaCB	138	80
2,2',3,4',5,6-HexaCB	147	80
2,2',3,4',5',6-HexaCB	149	80
2,2',4,4',5,5'-HexaCB	153	80
2,2',4,4',6,6'-HexaCB	155	80
2,3,3',4,4',5-HexaCB	156	80
2,3,4,4',5,6-HexaCB	166	80
3,3',4,4',5,5'-HexaCB	169	80
2,2',3,3',4',5,6-HeptaCB	177	80
2,2',3,4,4',5,5'-HeptaCB	180	80
2,2',3,4',5,5',6-HeptaCB	187	80
2,2',3,4',5,6,6'-HeptaCB	188	80
2,3,3',4,4',5,5'-HeptaCB	189	80
2,2',3,3',4,5,5',6'-OctaCB	199	80
2,2',3,3',5,5',6,6'-OctaCB	202	80
2,3,3',4,4',5,5',6-OctaCB	205	80
2,2',3,3',4,4',5,5',6-NonaCB	206	80
2,2',3,3',4,5,5',6,6'-NonaCB	208	80
DecaCB	209	80