

EO-5403

RoHS PBDE Cleanup Spike ($^{13}\text{C}_{12}$, 99%)1.2 mL
in Nonane

Labeled BDE	BDE#	Concentration (ng/mL)
4-MonoBDE ($^{13}\text{C}_{12}$, 99%)	3L	100
4,4'-DiBDE ($^{13}\text{C}_{12}$, 99%)	15L	100
2,4,4'-TriBDE ($^{13}\text{C}_{12}$, 99%)	28L	100
2,2',4,4'-TetraBDE ($^{13}\text{C}_{12}$, 99%)	47L	100
2,2',4,4',5-PentaBDE ($^{13}\text{C}_{12}$, 99%)	99L	100
2,2',4,4',5,5'-HexaBDE ($^{13}\text{C}_{12}$, 99%)	153L	200
2,2',4,4',5,6'-HexaBDE ($^{13}\text{C}_{12}$, 99%)	154L	200
2,2',3,4,4',5,6-HeptaBDE ($^{13}\text{C}_{12}$, 99%)	183L	200
2,2',3,4,4',5,6,6'-OctaBDE ($^{13}\text{C}_{12}$, 99%)	204L	200
2,2',3,3',4,4',5,6,6'-NonaBDE ($^{13}\text{C}_{12}$, 99%)	207L	500
DecaBDE ($^{13}\text{C}_{12}$, 99%)	209L	500