



Cambridge Isotope Laboratories, Inc.

3 Highwood Drive, Tewksbury, MA 01876 USA
PH: 800.322.1174 (N. America) PH: 978-749-8000
Fax: 978-749-2768 WEB: www.isotope.com

ES-5570

PFOS/PFOA Calibration Series CS1-CS5

5 x 0.25 mL in methanol

*All concentrations are in ng/mL***

Unlabeled	CS0.25*	CS1	CS2	CS3	CS4	CS5
Perfluorooctanoic acid (PFOA)	0.5	2	10	50	200	1000
Perfluorooctanesulfonate (PFOS), sodium salt	0.5	2	10	50	200	1000
Extraction						
Perfluorooctanoic acid (PFOA) (¹³ C ₈ , 99%)	50	50	50	50	50	50
Perfluorooctanesulfonate (PFOS), sodium salt (¹³ C ₈ , 99%)	50	50	50	50	50	50
Injection						
Perfluorononanoic acid (PFNA) (¹³ C ₉ , 99%)	50	50	50	50	50	50

**Not included in ES-5570 – available for separate purchase*

****Gravimetric concentration is based on the salt form, where applicable.**