



Cambridge Isotope Laboratories, Inc.

3 Highwood Drive, Tewksbury, MA 01876 USA

PH: 800.322.1174 (N. America) PH: 978-749-8000

Fax: 978-749-2768 WEB: www.isotope.com

ES-5653      Method 521 Nitrosamine Analyte PDS Mixture      1.2 mL in methylene chloride

<b>Unlabeled</b>	<b>Concentration (<math>\mu\text{g/mL}</math>)</b>
<i>N</i> -Nitrosodimethylamine	100
<i>N</i> -Nitrosomethylethylamine	100
<i>N</i> -Nitrosodiethylamine	100
<i>N</i> -Nitroso-di- <i>n</i> -propylamine	100
<i>N</i> -Nitroso-di- <i>n</i> -butylamine	100
<i>N</i> -Nitrosopyrrolidine	100
<i>N</i> -Nitrosopiperidine	100