

EC-4978

Method 1668A Labeled Clean-up Standard Solution

1.2 mL
in Nonane

Labeled PCBs	IUPAC*	Concentration ($\mu\text{g/mL}$)
2,4,4'-TriCB ($^{13}\text{C}_{12}$,99%)	28L	1.0
2,3,3',5,5'-PentaCB ($^{13}\text{C}_{12}$,99%)	111L	1.0
2,2',3,3',5,5',6-HeptaCB ($^{13}\text{C}_{12}$,99%)	178L	1.0