

EC-5370

EN-1948-4 PCB Sampling Standard ( $^{13}\text{C}_{12}$ , 99%)1.2 mL  
in Nonane

$^{13}\text{C}$ -Labeled Compounds	IUPAC#	Concentration (ng/mL)
$^{13}\text{C}_{12}$ -2,3,4,4'-TetraCB	60	100
$^{13}\text{C}_{12}$ -3,3',4,5,5'-PentaCB	127	100
$^{13}\text{C}_{12}$ -2,3,3',4,5,5'-HexaCB	159	100