

EC-4972

JIS PCB Type 2 Clean-up Standard Solution

1.2 mL  
in Nonane

<b>Labeled PCB</b>	<b>IUPAC*</b>	<b>EC-4972 Concentration (ng/mL)</b>	<b>EC-4972-A Concentration (µg/mL)</b>
3,3',4,4'-TetraCB ( <sup>13</sup> C <sub>12</sub> ,99%)	77	10	2
3,4,4',5-TetraCB ( <sup>13</sup> C <sub>12</sub> ,99%)	81	10	2
2,3',4,4',5-PentaCB ( <sup>13</sup> C <sub>12</sub> ,99%)	118	10	2
3,3',4,4',5-PentaCB ( <sup>13</sup> C <sub>12</sub> ,99%)	126	10	2
2,3,3',4,4',5-HexaCB ( <sup>13</sup> C <sub>12</sub> ,99%)	156	10	2
3,3',4,4',5,5'-HexaCB ( <sup>13</sup> C <sub>12</sub> ,99%)	169	10	2
2,3,3',4,4',5,5'-HeptaCB ( <sup>13</sup> C <sub>12</sub> ,99%)	189	10	2