

ED-935-A

Modified ED-935 TCDD Column Performance Check Solution  
(<sup>13</sup>C<sub>12</sub>, 99%; <sup>37</sup>Cl<sub>4</sub>, 96%)0.5 mL  
in Nonane

Labeled Component	Concentration (ng/mL)
2,3,7,8-Tetrachlorodibenzo-p-dioxin	100
1,2,3,4-Tetrachlorodibenzo-p-dioxin	100
1,4,7,8-Tetrachlorodibenzo-p-dioxin	100
1,2,3,7/1,2,3,8-Tetrachlorodibenzo-p-dioxin pair	100
1,2,7,8-Tetrachlorodibenzo-p-dioxin	200
<sup>13</sup> C <sub>12</sub> -2,3,7,8-Tetrachlorodibenzo-p-dioxin	250
<sup>37</sup> Cl <sub>4</sub> -2,3,7,8-Tetrachlorodibenzo-p-dioxin	7
<sup>13</sup> C <sub>12</sub> -1,2,3,4-Tetrachlorodibenzo-p-dioxin	500