

All concentrations are in ng/mL

Native BDEs	Congener**	CS1	CS2	CS3	CS4	CS5
2,4,4'-TriBDE	28	1	5	50	500	2500
2,2',4,4'-TetraBDE	47	1	5	50	500	2500
2,2',4,4',5-PentaBDE	99	1	5	50	500	2500
2,2',4,4',6-PentaBDE	100	1	5	50	500	2500
2,2',4,4',5,5'-HexaBDE	153	1	5	50	500	2500
2,2',4,4',5,6'-HexaBDE	154	1	5	50	500	2500
2,2',3,4,4',5,6-HeptaBDE	183	1	5	50	500	2500
DecaBDE	209	10	50	500	5000	25000
Labeled BDEs	Congener**	CS1	CS2	CS3	CS4	CS5
2,4,4'-TriBDE (¹³ C ₁₂ ,99%)	28L	100	100	100	100	100
2,2',4,4'-TetraBDE (¹³ C ₁₂ ,99%)	47L	100	100	100	100	100
2,2',4,4',5-PentaBDE (¹³ C ₁₂ ,99%)	99L	100	100	100	100	100
2,2',4,4',6-PentaBDE (¹³ C ₁₂ ,99%)	100L	100	100	100	100	100
2,2',4,4',5,5'-HexaBDE (¹³ C ₁₂ ,99%)	153L	100	100	100	100	100
2,2',4,4',5,6'-HexaBDE (¹³ C ₁₂ ,99%)	154L	100	100	100	100	100
2,2',3,4,4',5,6-HeptaBDE (¹³ C ₁₂ ,99%)	183L	100	100	100	100	100
DecaBDE (¹³ C ₁₂ ,99%)	209L	1000	1000	1000	1000	1000
Labeled Clean-up	Congener**	CS1	CS2	CS3	CS4	CS5
2,2',3,4,4',6-HexaBDE (¹³ C ₁₂ ,99%)	139L	100	100	100	100	100
Labeled Injection Internal	Congener**	CS1	CS2	CS3	CS4	CS5
2,2',5,5'-TetraCB (¹³ C ₁₂ ,99%)	PCB-52L	100	100	100	100	100
2,2',3,4,4',5'-HexaCB (¹³ C ₁₂ ,99%)	PCB-138L	100	100	100	100	100