

X

EC-4938

PCB Mixture-A

3 mL
in Nonane

Labeled PCBs	IUPAC*	Concentration ($\mu\text{g/mL}$)
3,3',4,4'-TetraCB ($^{13}\text{C}_{12}$,99%)	77	1.0
3,4,4',5-TetraCB ($^{13}\text{C}_{12}$,99%)	81	1.0
2',3,4,4',5-PentaCB ($^{13}\text{C}_{12}$,99%)	123	1.0
3,3',4,4',5-PentaCB ($^{13}\text{C}_{12}$,99%)	126	1.0
3,3',4,4',5,5'-HexaCB ($^{13}\text{C}_{12}$,99%)	169	1.0
2,2',3,4,4',5,5'-HeptaCB ($^{13}\text{C}_{12}$,99%)	180	1.0