

ECN-5490

PCN Cleanup Solution ($^{13}\text{C}_{10}$, 99%)5 mL
in Methanol:Isooctane

Labeled Component	PCN	Concentration (ng/mL)
1,2,3,4-TetraCN ($^{13}\text{C}_{10}$, 99%)	27	0.5
1,2,3,5,7-PentaCN ($^{13}\text{C}_{10}$, 99%)	52	0.5
1,2,3,4,5,7-HexaCN ($^{13}\text{C}_{10}$, 99%)	64	0.5
1,2,3,5,6,7-HexaCN ($^{13}\text{C}_{10}$, 99%)	67	0.5
1,2,3,4,5,6,7-HeptaCN ($^{13}\text{C}_{10}$, 99%)	73	0.5
OctaCN ($^{13}\text{C}_{10}$, 99%)	75	0.5