

ES-2044

PAH Surrogate Cocktail

1 mL  
50% Methylene chloride-D<sub>2</sub>/  
50% Methanol-OD

Labeled Component	Concentration (ug/mL)
Acenaphthylene (D <sub>8</sub> , 98%)	200
Benzo[ <i>a</i> ]pyrene (D <sub>12</sub> , 98%)	200
Benzo[ <i>g,h,i</i> ]perylene (D <sub>12</sub> , 98%)	200
Fluoranthene (D <sub>10</sub> , 98%)	200
Naphthalene (D <sub>8</sub> , 98%)	200
Phenanthrene (D <sub>10</sub> , 98%)	200
Pyrene (D <sub>10</sub> , 98%)	200