

ES-5321-5X10 Multi-Analyte Recovery Spiking Standard (¹³C, 99%)

5 X 10 mL
in 88% Hexane/
2% dodecane/
10% nonane

¹³ C-Labeled Component	IUPAC	Concentration (ng/mL)
1,2,3,4-TCDD (¹³ C ₆ , 99%)		2.5
2,2',3,3',4,5,5',6,6'-NonaCB (¹³ C ₁₂ , 99%)	208	10
3,3',4,4'-TetraBDE (¹³ C ₁₂ , 99%)	77	7.5
2,2',3,4,4',6-HexaBDE (¹³ C ₁₂ , 99%)	139	7.5