

ES-5474

CDC PCB Recovery standard for OH-PAHS ( $^{13}\text{C}_{12}$ , 99%)1.0 mL  
in Toluene

$^{13}\text{C}$ -labeled Component	IUPAC	Concentration (ng/mL)
2,4,4'-TriCB ( $^{13}\text{C}_{12}$ , 99%)	28	200
2,3,3',4,4'-PentaCB ( $^{13}\text{C}_{12}$ , 99%)	105	200
2,2',4,4',5,5'-HexaCB ( $^{13}\text{C}_{12}$ , 99%)	153	200
2,3',4,4',5,5'-HexaCB ( $^{13}\text{C}_{12}$ , 99%)	167	200